



## What are the benefits of UVR-M & UVR-Mi Air Purifiers?

The Biosan UV Air Purifiers are ideal for sanitising air in most indoor environments by efficiently destroying disease-causing agents (e.g., Viruses, Micro-organisms) by UV radiation.

- Full user protection from direct UV radiation
- Very low noise levels
- Simple to operate and maintain
- Digital control unit can track UV lamp operating time (UVR-Mi only)
- Easy to install – fix to a wall, or mount on the movable stand.



**\$1,159**  
UVR-Mi

**\$666**  
UVR-M

| SPECIFICATIONS                                  | UVR-M                         | UVR-Mi                        |
|---|-------------------------------|-------------------------------|
| UV radiation source: 25W Bactericidal (Philips) | 1 lamp                        | 2 lamps                       |
| UV radiation level                              | 18 mW / cm <sup>2</sup> / sec | 36 mW / cm <sup>2</sup> / sec |
| Air-flow productivity                           | 21m <sup>3</sup> / hour       | 29m <sup>3</sup> / hour       |
| Full user protection from direct UV light       |                               | Yes                           |
| Display   | —                             | LCD                           |
| UV lamp operation indicator                     |                               | Yes                           |
| UV lamp lifetime counter                        | No                            | Yes                           |
| Timer   | —                             | 1 min-24 hrs / non-stop       |
| Automatic switch ON/OFF                         | —                             | Yes                           |
| Lamp fault detection                            | —                             | Yes                           |

### ORDERING

| CATALOG №     | DESCRIPTION   | QTY. | PRICE (ex-GST)*   |
|---------------|---|------|-------------------|
| BS-040110-AA3 | Biosan UVR-Mi Air Purifier / UV Cleaner-Recirculator, (2 lamps), 110 x 135 x 660 mm | 1    | <b>\$1,159.00</b> |
| BS-040105-AA3 | Biosan UVR-M Air Purifier / UV Cleaner-Recirculator, (1 lamp), 110 x 135 x 660 mm   | 1    | <b>\$666.00</b>   |
| BS-040105-AK  | Biosan UVR-S tripod stand, 510 x 510 x 1100 mm                                      | 1    | <b>\$222.00</b>   |

**ORDER NOW**

Email [info@fisherbiotec.com](mailto:info@fisherbiotec.com) for more details or visit our [website](#)